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The Shin-Etsu products “SILICOSEN” and “BIO-SILICO” are culture plugs specially processed in silicone rubber with continuous bubbles. The products have excellent properties such as permeability, operability, heat resistance and chemical resistance.

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Shin-Etsu SILICOSEN

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- 1.** 綿栓に比べ取り扱いが容易で、省力化、人件費の大幅削減が可能です。  
Easy to use in comparison to cotton plugs which means that massive savings can be made in labor and labor costs.
- 2.** 綿栓に比べ水分蒸発速度が1/2以下です。培養液の異常濃縮、PH変化が少なくなります。  
Moisture evaporation rate is less than half of cotton plugs which means that heteropycnosis of the culture solution and PH changes are reduced.
- 3.** 耐熱性に優れ乾熱滅菌(180℃)、高圧蒸気滅菌による繰り返し使用が可能です。  
Excellent heat resistance property which means that plugs can be dry sterilized (180℃) and high-pressure steam sterilized for repeated use.
- 4.** 耐薬品性に優れています。  
Excellent chemical resistance property.
- 5.** 熱湯洗浄、消毒が可能です。また中性洗剤で洗浄できます。  
Can be hot-water washed and disinfected. Also can be washed in neutral detergent.
- 6.** 疎水性のため、容器とのシール部分での菌の繁殖を抑制します。  
Hydrophobic property means that bacteria proliferation is suppressed at the sealing area with the container.

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Amount of moisture evaporation is minimal, which means that there is less drying of the culture medium.

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Possesses superb permeability and is type where importance is placed on mounting quality. Suitable for shaking cultures.

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Possesses superb permeability and suitable for aerobic bacterium cultures and shaking cultures.

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





Moisture evaporation rate is half or less of cotton plug for T type and approximately the same as cotton plug for C type.

Test method:

Inject 5ml of water into a medium-sized test tube, seal tube with culture plug and leave standing. The fluctuation in weight due to moisture evaporation is measured with the elapse of time.

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SILICOSEN and BIO-SILICO prevents the stray entry of airborne bacteria for both stationary cultures and shaking cultures.

\*Number of test tubes with stray entries of airborne bacteria / number of submitted test tubes.

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#### Comparison of Features among Culture Plug Types

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A. Optimum, B. Suitable, C. Not suitable

※ Varies depending on style of cotton plug manufacture and usage.

	T-タイプ T-Type	S-タイプ S-Type	バイオシリコ BIO-SILICO	綿 栓 Cotton Plug	紙 栓 Paper Plug
栓の着脱が容易 Easy fitting	A	A	A	C※	A
植菌操作性 Inoculation	B	B	B	B	B
繰り返し使用 Repeated Use	B	B	B	C	C
洗浄・乾燥 Easy washing, Quick dryup	B	B	A	C	C
容器への保持 Fitting to tubes, flasks	A	A	A	C※	A
高圧蒸気滅菌 Autoclaving	B	B	B	B	B
乾熱滅菌 Dry air sterilization	A	A	B	B	B
火焰滅菌 Flame sterilization	B	B	C	A	A

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×: >75%

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Test method: Immerse and leave for one week in various types of chemicals at room temperature (23℃) and then observe (evaluate) tensile strength changes.

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Proper Sterilization Methods for SILICOSEN and BIO-SILICO

	シリコセン SILICOSEN	バイオシリコ BIO-SILICO
高圧蒸気滅菌 (121°C 20min) High-pressure steam sterilization	可能 (最適滅菌法) Possible (optimum sterilization method)	可能 (最適滅菌法) Possible (optimum sterilization method)
乾熱滅菌① (120°C 20~30min) Dry air sterilization	可能 Possible	可能 Possible
乾熱滅菌② (180°C 15~30min) Dry air sterilization	可能 Possible	不可 Not possible
火焰滅菌 Flame sterilization	瞬時可能 Possible for an instant	不可 Not possible

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Items of caution concerning high-pressure steam sterilization

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Proper Washing Methods for SILICOSEN and BIO-SILICO

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Methods include rough washing in cold or hot water followed by boiling or ultrasound washing. Also neutral detergent is the suitable detergent. But note that all detergent must be completely rinsed off.

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In Cases where Surface Culture Medium Hardens and Blackens with Mold

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If mold has infiltrated, carbonizes and blackens the plug so that even washing in hot water with neutral detergent does not remove the mold, replace plug with a new one.

Note that plugs are to be disposed of as plastic and, as such, should be handled accordingly.

- These products are ventilated plugs for use in cultivating cultures and growths. Do not use them for any other application.
- Always confirm the safety of intended use by checking the Shin-Etsu application details prior to using the products.
- Select product size that fits mouth diameter of container. Forcibly insertion of a SILICOSEN or BIO-SILICO plug that is too large into a container causes loss of ventilation and leads to trouble.
- The products are suitable for some cultures and not for others. We recommend that a test culture run be performed prior to actual use.

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- The data and information presented in this catalog may not be relied upon to represent standard values. Shin-Etsu Polymer reserves the right to change such data and information, in whole or in part, in this catalog, including product performance standards and specifications without notice.
- Users are solely responsible for making preliminary tests to determine the suitability of products for their intended use. Statements concerning possible or suggested uses made herein may not be relied upon, or be construed, as a guaranty of no patent infringement.
- Our silicone products described herein have been designed, manufactured and developed solely for general industrial use only; such silicone products are not designed for, intended for use as, or suitable for, medical, surgical or other particular proposes. Users have the sole responsibility and obligation to determine the suitability of the silicone products described herein for any application, to make preliminary tests, and to confirm the safety of such products for their use.
- Users must never use the silicone products described herein for the purpose of implantation into the human body and/or injection into humans.
- Users are solely responsible for exporting or importing the silicone products described herein, and complying with all applicable laws, regulations, and rules relating to the use of such products. Shin-Etsu Polymer recommends checking each pertinent country's laws, regulations, and rules in advance, when exporting or importing, and before using the products.
- Please contact Shin-Etsu Polymer Group III Sales&Marketing Division II before reproducing any part of this catalog.

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The Development and Manufacture of silicone products are based on the following registered international quality, environmental, medical devices - quality management standards.

Kodama Plant	ISO9001 (2681603)	ISO14001 2691359	ISO13485 DNKFRC95078A)
Nagano Plant	ISO9001 (3014595)	ISO14001 JCQA-E-0056	ISO13485 DNKFRC96822A)
Shin-Etsu Polymer (Malaysia) Sdn. Bhd.	ISO9001 (Q209087)	ISO14001 E209041	ISO13485 BVC130000)

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Shin-Etsu Polymer Co., Ltd.
**Shin-Etsu Polymer Singapore Pte. Ltd.**
**Shin-Etsu Polymer Hong Kong Co., Ltd.**